



More information on the website  
[radwag.com/en/info,w1,9XE](http://radwag.com/en/info,w1,9XE)

# 100 $\mu$ l $\div$ 1000 $\mu$ l - RW8-106-20-9 Automatic Pipettes



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
Volume	100 $\mu$ l $\div$ 1000 $\mu$ l
Systematic error	2.0% / 0.6%
Random error	0.7% / 0.2%
Tips	1000 $\mu$ l



## Compatible with

### XA 5Y.M.A.P Microbalance

- XA 6/21.5Y.M.A.P Microbalance [WL-112-1000]
- XA 21.5Y.M.A.P Microbalance [WL-112-0001]
- XA 53.5Y.M.A.P Microbalance [WL-112-0003]
- XA 21/52.5Y.M.A.P Microbalance [WL-112-1001]
- XA 52.5Y.M.A.P Microbalance [WL-112-0002]

### Microbalance for Pipette Calibration MYA 5Y.P

- MYA 21.5Y.P Microbalance [WL-109-0023]

### Automatic Device for Multichannel Pipette Calibration

- AP-12.5Y Automatic Device for Multichannel Pipette Calibration [WL-101-0416]
- AP-12.1.5Y Automatic Device for Multichannel Pipette Calibration [WL-101-0423]