



More information on the website
radwag.com/en/info,w1,UAW

0.5 µl ÷ 10 µl - RW8-101-20-9 Automatic Pipettes



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
Volume	0.5 µl ÷ 10 µl
Systematic error	2.5% / 1.0%
Random error	1.5% / 0.8%
Tips	10 µl



Compatible with

XA 5Y.M.A.P Microbalance

- XA 6/21.5Y.M.A.P Microbalance [WL-112-1000]
- XA 21.5Y.M.A.P Microbalance [WL-112-0001]
- XA 53.5Y.M.A.P Microbalance [WL-112-0003]
- XA 21/52.5Y.M.A.P Microbalance [WL-112-1001]
- XA 52.5Y.M.A.P Microbalance [WL-112-0002]

Microbalance for Pipette Calibration MYA 5Y.P

- MYA 21.5Y.P Microbalance [WL-109-0023]

Automatic Device for Multichannel Pipette Calibration

- AP-12.5Y Automatic Device for Multichannel Pipette Calibration [WL-101-0416]
- AP-12.1.5Y Automatic Device for Multichannel Pipette Calibration [WL-101-0423]